

BECAUSE
TIME
MATTERS

IDYLLA™ EGFR MUTATION ASSAY



Idylla™ EGFR Mutation Assay is for research use only. Not for use in diagnostic procedures.

THE IDYLLA™ EGFR MUTATION ASSAY

FAST & RELIABLE INFORMATION ON EGFR MUTATION STATUS



Comprehensive assay covering 51 mutations in EGFR exons 18, 19, 20 and 21



Directly from 1 FFPE tissue section



150 minutes assay turn around time



Highly sensitive & standardized platform



< 2 minutes hands-on time




Fully automated solution - suitable for any lab

IDYLLA™ EGFR MUTATION ASSAY TISSUE SPECIMEN REQUIREMENTS

One of the biggest challenges in oncology biomarker testing is the ability to obtain samples of sufficient size and quality. The Idylla™ assay is designed to work with a minimal amount of sample:

- 1 x 5 µm FFPE tissue section
- Neoplastic cell content ≥ 10%
- If < 10%, macrodissection needed



IDYLLA™ EGFR MUTATION ASSAY EXCELLENT PERFORMANCE



The Idylla™ EGFR Mutation Assay shows high concordance to reference methods across a variety of different sample types

| STUDY | # SAMPLES | REFERENCE METHOD | SAMPLE TYPE | CONCORDANCE |
|--|-----------|------------------|--|-------------|
| Brohawn et al. 2018 ¹ | 23 | PCR | FFPE tissue | 100% |
| De Luca et al. 2018 ² | 43 | Ion Torrent NGS | DNA | 100% |
| De Montpréville et al. 2017 ³ | 93 | Sentosa NGS | FFPE tissue, fresh frozen tissue, cytological samples | 100% |
| Arcila et al. 2018 ⁴ | 62 | NGS | FFPE tissue, DNA, cytological samples (cell pellets, smears), touch preps, NGS libraries | 100% |
| Al-Turkmani et al. 2018 ⁵ | 8 | Ion Torrent NGS | Fresh touch preps | 100% |
| Al-Turkmani et al. 2018 ⁶ | 34 | Ion Torrent NGS | FFPE tissue | 100% |



The Idylla™ EGFR Mutation Assay can significantly reduce the time until the biomarker test result becomes available

| STUDY | REFERENCE METHOD | TIME TO RESULT VS. REFERENCE METHOD |
|------------------------------------|------------------|--|
| Brohawn et al. 2018 ¹ | In-house PCR | < 160 minutes vs. 35 hours 5 minutes (13 times faster than reference method) |
| Ghigna et al. 2017 ⁷ | Send-out NGS | ≤ 1 week (88%) / same day (37%) vs. ≤ 3 weeks (69%) / > 8 days (100%) |
| Mackinnon et al. 2019 ⁸ | In-house NGS | 1-2 days vs. 2-15 days* |

* Data shown as average for Idylla™ KRAS, BRAF, NRAS-BRAF and EGFR Mutation Assays



The Idylla™ EGFR Mutation Assay is able to analyze samples which do not have sufficient quantity or quality for NGS testing or have failed on previous NGS testing platforms

| STUDY | REFERENCE METHOD | VALID RESULTS IN % (# OF SAMPLES) VS. REFERENCE METHOD |
|--------------------------------------|------------------|---|
| De Luca et al. 2018 ² | Ion Torrent NGS | 93% (63/68) vs. 63% (43/68) |
| Al-Turkmani et al. 2018 ⁶ | Ion Torrent NGS | 98% (39/40) vs. 85% (34/40) |
| Mackinnon et al. 2019 ⁸ | Ion Torrent NGS | 25 NGS QNS samples: valid Idylla™ results in all cases, including actionable results in 8 cases |

IDYLLA™ EGFR MUTATION DETECTION



EGFR is a key component regulating tumor cell proliferation and growth and it is frequently mutated in different types of human cancers including lung cancer.⁹

The Idylla™ EGFR Mutation Assay provides qualitative detection of exon 18 (G719A/C/S), exon 21 (L858R, L861Q), exon 20 (T790M, S768I) mutations, exon 19 deletions and exon 20 insertions in the EGFR oncogene.



| | | |
|----------------|----------|--|
| Exon 18 | G719A | c.2156G>C |
| | G719S | c.2155G>A |
| | G719C | c.2155G>T; c.2154_2155delinsTT |
| Exon 19 | Del 9 | c.2238_2248delinsGC; c.2239_2248delinsC; c.2240_2248del; c.2239_2247del |
| | Del 12 | c.2239_2251delinsC; c.2240_2251del |
| | Del 15 | c.2235_2249del; c.2236_2250del; c.2239_2253del; c.2240_2254del; c.2238_2252del; c.2237_2251del; c.2235_2252delinsAAT; c.2237_2252delinsT; c.2234_2248del; c.2236_2253delinsCTA; c.2237_2253delinsTA; c.2235_2251delinsAG; c.2236_2253delinsCAA; c.2230_2249delinsGTCAA |
| | Del 18 | c.2240_2257del; c.2237_2255delinsT; c.2239_2256del; c.2236_2253del; c.2239_2258delinsCA; c.2237_2254del; c.2238_2255del; c.2237_2257delinsTCT; c.2236_2255delinsAT; c.2236_2256delinsATC; c.2237_2256delinsTT; c.2237_2256delinsTC; c.2235_2255delinsGGT |
| | Del 21 | c.2238_2258del; c.2236_2256del |
| | Del 24 | c.2253_2276del |
| Exon 20 | T790M | c.2369C>T |
| | S768I | c.2303G>T |
| | insG | c.2310_2311insGGT |
| | insASV9 | c.2307_2308insGCCAGCGTG |
| | insASV11 | c.2309_2310delinsCCAGCGTGGAT |
| Exon 21 | insSVD | c.2311_2312insGCGTGGACA |
| | insH | c.2319_2320insCAC |
| | L858R | c.2573T>G; c.2573_2574delinsGT; c.2573_2574delinsGA |
| | L861Q | c.2582T>A |

THE IDYLLA™ ADVANTAGE



The fully automated Idylla™ EGFR Mutation Assay provides fast and reliable information on EGFR mutation status directly from FFPE tissue.¹⁻⁷



The Idylla™ EGFR Mutation Assay shows high concordance to reference methods across a variety of different samples types.¹⁻⁶



The Idylla™ EGFR Mutation Assay provides high accuracy with low sample input showing a reduced failure rate compared to reference technologies.^{2,6,8}

The fast turnaround time and ease of use makes the Idylla™ solution suitable for any lab - as a cost effective, quick and accurate in-house system for small centers that lack highly trained staff and molecular expertise as well as a reliable STAT solution complementing comprehensive NGS profiling at larger centers.



REFERENCES

- (1) Brohawn DG et al. (2019) Journal of Diagnostic Techniques and Biomedical Analysis; 7:2.
- (2) De Luca et al. (2018) Journal of Clinical Pathology; 71:745-750.
- (3) De Montpréville et al. (2017) Pathology - Research and Practice; 213,7:793-798.
- (4) Arcila M et al. (2018) Journal of Molecular Diagnostics; 20:6. Abstract #ST027.
- (5) Al-Turkmani M (2018) Journal of Molecular Diagnostics; 20:6. Abstract # TT060.
- (6) Al-Turkmani M (2018) Journal of Molecular Diagnostics; 20:6. Abstract # ST125.
- (7) Ghigna M et al. (2018) Journal of Thoracic Disease; 10(7):4653-4658.
- (8) Mackinnon A et al. (2019) Accepted for publication in Journal of Molecular Diagnostics.
- (9) Sigismund S et al. (2018) Molecular Oncology; 12(1):3-20.

Contact details

Biocartis US Inc.
2 Pierce Place Suite 1510
Itasca, IL 60143
1-844-4-IDYLLA (1-844-443-9552)
www.biocartis.com/US
customerserviceUS@biocartis.com

Follow us on   

For Research Use Only. Not for use in diagnostic procedures.

Copyright information

The Idylla™ platform is available as an IVD in Europe and many other countries outside the US and is cleared in the US under K163628. Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis. © January 2023, Biocartis US Inc. All rights reserved.