

IDYLLA™ CONNECTIVITY AND SOFTWARE APPLICATIONS

Extend your user experience by integrating the Idylla™ Platform in your laboratory network



Access and print the Idylla™ Test Result Report from any computer in your lab



Get automatic access to the latest software updates in one click



Upload the Idylla™ Test Results into your LIS system (HL7 compliant)



Gain insights into your run data with Idylla™ Explore or Idylla™ Visualizer



Receive immediate remote troubleshooting and service



MORE INSIGHTS INTO YOUR MOLECULAR TEST DATA WITH IDYLLA™ EXPLORE OR IDYLLA™ VISUALIZER

Access Idylla™ Test Results as presented on the report

Filter Idylla™ Target Results to display all valid and invalid PCR curves

Generate a PDF curve report

Visualize the PCR curves from Idylla™ Test Results

Access Cq and ΔCq values for each target

HEAD-TO-HEAD COMPARISON

	Idylla™ Explore	Idylla™ Visualizer
Access	Any computer around the globe	Local computer
Computer Specifications	Web browser	Windows 10 or higher - 64 bit 10 GB or more available free disk space 4 GB or more RAM 1920x1080 or higher display resolution
Run data export	Microsoft Excel 2013 or higher	Microsoft Excel 2013 or higher
Report viewer	PDF viewer	PDF viewer
Installation	N/A	Local
Console data transfer	Idylla™ Connect	USB, Direct RJ45 Ethernet cable or Hospital LAN
Supported Tests	IVD* and RUO	IVD* and RUO
Data storage	Cloud	Computer
Upgrade by	Biocartis	Customer
Computer by	Customer	Customer
Compliance	All regulatory frameworks and associated industry standards are described in the Idylla™ Data Security & Compliance Guide.	

Contact your local Biocartis representative to receive a proposal that is tailored to your laboratory set-up and needs.

Biocartis NV
 Generaal De Wittelaan 11B
 2800 Mechelen - Belgium
 +32 15 632 888

Biocartis US Inc.
 2 Pierce Place Suite 1510
 Itasca, IL 60143 - USA
 1-844-4-IDYLLA (1-844-443-9552)

Follow us on www.biocartis.com
customerservice@biocartis.com
customerserviceUS@biocartis.com

* In markets where available.
 Idylla™ Platform is CE-marked in Europe in compliance with EU IVD Regulation 2017/746, cleared for sales in the US and registered in many others countries. Idylla™ Explore is for Research Use Only. Not for use in diagnostic procedures. Idylla™ is available for sale in Europe, the US and many other countries. Please check availability with a Biocartis representative. Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis NV. © October 2023, Biocartis NV. All rights reserved.

