

BECAUSE
TIME
MATTERS

IDYLLA™ ctEGFR MUTATION ASSAY



Idylla™ ctEGFR Mutation Assay is for research use only. Not for use in diagnostic procedures.

FAST AND EASY ctEGFR TESTING



COMPREHENSIVE ASSAY

- Covers 49 mutations including insertions and deletions in exon 18, 19, 20 and 21 in the *EGFR* gene



EASE OF USE

- Fully automated sample-to-result process
- Walk-away system
- All PCR reagents integrated in a single cartridge
- Contamination-controlled design



FAST RESULTS

- ± 2 minutes hands-on time
- Short turnaround time ± 160 minutes from sample to result



SUITABLE FOR ANY LAB

- No special expertise required
- Minimal infrastructure needed



LIQUID BIOPSY

- Directly from 2 ml of plasma



ACCESSIBLE

- Access on demand
- No need for pre-processing or batching



THE IDYLLA™ ctEGFR MUTATION ASSAY SHOWS HIGH CONCORDANCE TO NGS¹

The Idylla™ ctEGFR Mutation Assay* compared to Next Generation Sequencing, shows **high sensitivity** and **ease of use** combined with **fast turnaround time** for the testing of 49 relevant EGFR mutations in 2 ml of NSCLC plasma samples.

- Thirty-four mutations were detected by NGS and confirmed by the Idylla™ ctEGFR prototype assay.
- An additional 8 mutations were detected by the Idylla™ ctEGFR prototype assay and confirmed by Cobas.

MUTATIONS DETECTED

Idylla™
42

NGS
34



(1) Reijans et al., ESMO 2019
*Prototype


COMPLEMENT YOUR TISSUE TESTING

Obtaining samples for biomarker testing of sufficient quantity and quality is known to be a major challenge in several cancer types including lung cancer:

- Tissue samples might be too small
- Tumor content might be too low
- Tumor might be hard to reach

Especially for these cases, testing plasma samples can be a rapid and cost-effective alternative and a valuable complement to tissue testing.

LIQUID BIOPSY WORKFLOW IDYLLA™ ctEGFR MUTATION ASSAY



Contact details
Biocartis US Inc.
2 Pierce Place Suite 1510
Itasca, IL 60143
1-844-4-IDYLLA (1-844-443-9552)
www.biocartis.com/US
customerserviceUS@biocartis.com

Follow us on   

For Research Use Only. Not for use in diagnostic procedures.

Copyright information

The Idylla™ platform is available as an IVD in Europe and many other countries outside the US and is cleared in the US under K163628. Biocartis and Idylla™ are registered trademarks in Europe, the US and many other countries. The Biocartis and Idylla™ trademarks and logos are used trademarks owned by Biocartis. © January 2023, Biocartis US Inc. All rights reserved.